

## Intensity Normalization

# Intensity normalization

- ▶ Conventional MRI intensities (T1-w, T2-w, PD, FLAIR) are acquired in arbitrary units
- ▶ Images are not comparable across scanners, subjects, and visits, even when the same protocol is used.
  - ▶ This affects algorithm performance, prediction, inference.
  - ▶ Even simple things like thresholding an image
- ▶ Intensity normalization brings the intensities to a common scale across people.
- ▶ We will discuss subject-level intensity normalization.


# Goals of this tutorial

- ▶ Visualize intensity distributions from different subjects and tissue classes.
  - ▶ Because MS lesions are in the white matter, comparable white matter distributions are the target goal
- ▶ Apply the White Stripe intensity normalization (Shinohara et al. 2014).




# Visualizing whole-brain intensities

- ▶ For the moment, we will work with the T1-w images from the training data.
- ▶ Full brain densities are mixtures of the three tissue class distributions.

**Distribution of all Voxels in Brain Mask**



# Visualizing the intensities by tissue class



## Whole-brain normalization

- ▶ Let's Z-score each voxel using mean  $\mu_{WB}$  and standard deviation  $\sigma_{WB}$  computed from all voxels in the brain mask.


$$T1^{WB} = \frac{T1 - \mu_{WB}}{\sigma_{WB}}$$

- ▶ `zscore_img` is a function in `neurobase` that does this.
- ▶ It takes an image and a binary mask. The default is to use all voxels in the brain mask.

```
zscore_img(img = img, mask = mask)
```

# Whole-brain normalized intensities

**CSF Before**




**CSF After**



## Whole-brain normalized intensities

- Gray matter distributions are more comparable.

### Gray Matter Before




### Gray Matter After



## Whole-brain normalized intensities

- ▶ White matter distributions are more comparable.

### White Matter Before



### White Matter After

## Other Normalizations: White Stripe

- ▶ Whole-brain normalization may be sensitive to outliers.
- ▶ Lesions in MS can have very high intensities, which lead to bad estimates of mean/variance
  - ▶ Other more robust transformations may be used, such as using the median to center, IQR to scale, etc.
- ▶ White Stripe (Shinohara et al. 2014) is based on parameters obtained from a sample of normal appearing white matter (NAWM), which is robust to outliers.
  - ▶ The idea is to make normal appearing white matter comparable across subjects and visits.

# White Stripe normalization

Procedure 1. Find white matter area on histogram

2. Estimate mean  $\mu_{WS}$  and variance  $\sigma_{WS}$  of voxel intensities in that area
3. Normalize with those means/variances:

$$T1^{WS} = \frac{T1 - \mu_{WS}}{\sigma_{WS}}$$

- ▶ After normalization, NAWM will have a standard normal distribution and units will be in standard deviations of NAWM.
- ▶ Gray matter and CSF distributions may not be comparable after White Stripe.

# White Stripe normalization




Figure 1: whitestripe


## White Stripe normalization code

```
ind = whitestripe(img = t1, type = "T1", stripped = TRUE)$v  
ws_t1 = whitestripe_norm(t1, indices = ind)
```

- ▶ The `whitestripe` function takes an image, image type (in our case T1), and a logical indicating whether the image has been skull stripped.
- ▶ The indices of voxels in the NAWM used for estimating the normalization parameters are located in the list element `$whitestripe.ind`.
- ▶ The function `whitestripe_norm` takes an image and the indices from a call to `whitestripe` and returns the White Stripe normalized image as a `nifti`.

# WhiteStripe normalized intensities

**Whole-brain: CSF**




**White Stripe: CSF**



# WhiteStripe normalized intensities

**Whole-brain: White Matter**



Whole-brain Normalized Intensities

**White Stripe: White Matter**



# Conclusions

- ▶ Intensity normalization is an important step in any image analysis with more than one subject or time point to ensure comparability across images.
- ▶ White Stripe normalization may work better and have better interpretation than whole-brain normalization for subsequent lesion segmentation algorithms and analysis.
- ▶ Other intensity normalization methods that make intensities comparable across subjects for all tissues exist.
  - ▶ RAVEL, which is an extension of WhiteStripe is one example (Fortin et al. 2016).
  - ▶ Located at <https://github.com/Jfortin1/RAVEL>
  - ▶ This was shown to have better comparability than histogram matching

# References

Fortin, Jean-Philippe, Elizabeth M Sweeney, John Muschelli, Ciprian M Crainiceanu, Russell T Shinohara, Alzheimer's Disease Neuroimaging Initiative, and others. 2016. "Removing Inter-Subject Technical Variability in Magnetic Resonance Imaging Studies." *NeuroImage* 132. Elsevier: 198–212.

Shinohara, Russell T, Elizabeth M Sweeney, Jeff Goldsmith, Navid Shiee, Farrah J Mateen, Peter A Calabresi, Samson Jarso, et al. 2014. "Statistical Normalization Techniques for Magnetic Resonance Imaging." *NeuroImage: Clinical* 6. Elsevier: 9–19.